

# Infu-Surg®

## Disposable Pressure Infuser Bag

Ethox International presents the Infu-Surg® disposable pressure infuser bag. The Infu-Surg® provides reliable results and quality performance, but unlike traditional means of IV infusion, the Infu-Surg® is equally convenient, beneficial, and simple to operate.

These infuser bags are individually packaged and designed for single-patient use, which essentially prevents complications during important hospital transport. The bags stay with a patient as they are moved into different departments, saving time and confusion. Disposable products encourage productivity as well. Time is saved for other crucial tasks because cleaning and sterilizing procedures can be eliminated. Most importantly, single-patient use of these devices prevents any possibility for cross contamination.

The bags load quickly and easily, and are single-hand operated, which allows the health care worker to simultaneously control other operating devices. The Infu-Surg® is adaptable with most IV bags and is available in multiple sizes. In addition, the Infu-Surg® has wide versatility because the same methods and principles can be applied to the infusion of a variety of fluids. With features such as a color-coded gauge, a pressure relief valve, and a single-hand stopcock - the Infu-Surg® provides superior reliability, cost-effectiveness, and safety.



### Product Features

- **Simple Use**
  - Single-handed inflation or deflation
  - Three-way stopcock ensures precise control for adjustments
  - No external power source is required
- **Completely Disposable**
  - Eliminates cleaning and repair
  - Prevents any cross contamination
- **Convenient**
  - Color-Coded Gauge for quick feedback on pressure status
  - Pressure Relief Valve protects the bag from over-inflation
  - Designed with a transparent front panel for observation
  - External inflation capabilities
- **Available in 500cc, 1000cc, and 3000cc Sizes**

### ORDER INFORMATION

Catalog #	Product	Quantity	Weight
4005	500cc Model	25 per case (or box of 5)	9 lbs
4010	1000cc Model	25 per case (or box of 5)	9 lbs
4030	3000cc Model	10 per case (or box of 5)	9 lbs